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INSTRUCTION FOR USE FOR ABUTMENTS BIO-1®

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1. **DESCRIPTION** : Healing abutment or standard / standard maxillary .

Abutment for dental implant manufactured, directly connected to the implant endo-osseous, designed to facilitate the prosthetic restoration.

2. **INDICATIONS** : Restore all or part of the masticatory function of a patient. The abutment makes it possible to create a dental prosthesis on a dental implant already in place. And allows to treat partial and total edentulousness. It is compatible with the following material :

Biosmile connection

Brånemark NR, RP and WP external hexagonal connection . Healing abutment diam 4.8 and 5.8. Abutment multiunit diam 4.8 and 6.0 with 0 °, 17 ° and 30 ° angle

Straumann Tissue Level RN and WN Connection

Straumann Bone Level NC and RC connection

Nobel CC conical connection

Nobel Replace trilob connection

Astra Tech OsseoSpeed and EV connection

Antogyr Axiom connection

Login Abutment Multi-Unit NR / PR diam diam WP 4.8 and 6.0

Zimmer Tapered Screw -Vent connection diam 3.5, 4.5 and 5.7

Connection Biotech dental diam 4.0 and 5.0

Global D In- Kone connection diam 4.0 and 5.0

There are 4 types of abutment :

1 ° Healing : as part of the implant placement, this abutment must be placed temporarily during the healing of the gum. Do not put false teeth on it, it goes in the maxillary part of the jaw as it does not interfere with the patient's aesthetic.

2 ° Temporary : same as the healing abutment but this is temporarily part incisor or lateral in order to place a false tooth during healing. This improves the aesthetics of the patient.

3 ° Standard : this is the final abutment, after healing.

4 ° Multiunit / rotational : they are used exclusively for the placement of a serie of implants to replace several contiguous teeth like a " bridge ".

To ensure the compatibility and the right choice , contact BIOSMILE INTEGRATION .

3. **CONTRAINDICATIONS** : The practitioner must perform a prior medical examination to ensure that the patient does not present any contraindications and in order to take into consideration the potential failure factors, a preoperative evaluation of the patient is necessary to determine any factor likely to expose the patient to a risk resulting from the use of a product. The products should not be used in patients considered medically unfit for conventional oral surgery. The products must not be used in patients for whom the imperatives of the number, position and dimensions of the implants cannot be respected in order to obtain a support adapted to the loads. Check if the patient has an allergy or hypersensitivity to titanium.

4. **PROTOCOL** : Any modification of the abutment must be done outside the oral cavity. Attention, due to the size of the abutments and screws, inform that the patient could inhale or ingest the product. To connect and to tighten the abutment and make sure there is no occlusal interference. It is recommended to carry out an X-ray check in order to ensure the correct positioning of the abutment. In order to achieve a harmonious and natural gum contour, a minimum distance of 1.25 mm will be required between the edge of the healing abutment and the adjacent teeth. The healing abutments should extend approximately 1 mm above the gumline at the highest point.

To tighten the abutment with the screwdriver and the Torque wrench, choosing the right torque :

- Nobel : 35Ncm
- Astra : 25Ncm for the final abutments, 15 for the healing abutments
- Straumann : 35Ncm
- Antogyr : 35Ncm
- Zimmer : 30Ncm
- Global D : 15 Ncm
- Biosmile : 35Ncm
- Multiunit : 15Ncm

Surgical protocols and dental operative techniques are widely documented by the original manufacturers of dental implants . If in doubt about the use of a product, contact BIOSMILE INTEGRATION ®.

5. **SIZE SELECTION** : There is the RP (regular platform) , NP (narrow platform) and the WP (wide platform) , chosen by the practitioner according to the implant placed and his size. If in doubt about the use of a abutment product , contact BIOSMILE INTEGRATION ®.

6. **CLEANING AND STERILIZATION** :

Before to use in the patient's mouth, non-sterile devices should be sterilized in a gravity steam autoclave. Only one device should be sealed in packaging for appropriate sterilization. BIO-1® products are compatible with approved sterilization methods commonly used in dental practice : 132-136 ° C for 3 minutes.

If in doubt about the sterilization process, contact BIOSMILE INTEGRATION ® . The abutment BIO-1® are disposable .

7. **STORAGE AND HANDLING** : An aseptic environment must be used for the use of the products.

8. **DESCRIPTION AND GENERAL PRECAUTIONS** : The products must be used by practitioners with a good knowledge of surgical and dental prosthetic techniques. Each patient should be carefully examined and assessed to determine their radiographic, psychological and physical condition. Close cooperation between the surgeon, the restorative dentist and the dental technician is essential for the success of the treatment. Practitioners are advised to follow theoretical and practical training in order to acquire the appropriate techniques, including biomechanical criteria and radiographic evaluation. Of improper techniques can cause failure of the prosthesis implant or implant itself.

9. **OPERATING PRECAUTIONS** : Every effort should be made to minimize damage to the host tissue, taking particular care to avoid thermal and surgical trauma and to eliminate all contaminants and sources of infection. Compared to conventional oral surgery, the surgical protocol requires great precision and attention, and the margins for acceptable tissue preparation are much narrower. All instruments used in surgery must be maintained in good condition. Care should also be taken to ensure that the instruments do not damage the components. Due to the small dimensions of the components, care should be taken to

prevent the patient from sucking or swallowing them. The user must prevent the risks of inhalation / ingestion of the product by the patient. After fitting the prosthesis, the surgeon will assess the quality and quality of the tissues and the initial stability of the prosthesis in order to determine the moment at which the products can be put into operation, In particular for the prosthetic parts, the correct distribution of the stresses is particularly important: passive adaptation and placement of the bridge on the implant abutments; adjustment of the occlusion relative to the opposing jaw; elimination of excessive transverse forces, especially in cases of immediate loading.

10. **INSTALLATION WARNING** : Each implant system has specific design features for the abutments and prosthetic components.


Combining components not configured or sized to accommodate these specificities can cause mechanical failure, tissue damage, or unsatisfactory aesthetic results.

If in doubt about the compatibility of BIO-1® products with other dental implant systems and the methods of using them, contact BIOSMILE INTEGRATION® directly . Preoperative soft and hard tissue deficiencies can compromise the aesthetic result.

Side effects: dental surgical techniques are accompanied by contraindications and risks of use. These are widely documented in the dental literature.

The BIO-1® implant system is composed by specific complementary parts that cannot be replaced by any other manufacture than that of BIOSMILE INTEGRATION® .

Any use of prosthetic material and / or screws not validated by BIOSMILE INTEGRATION ® entails, de facto, the forfeiture of the manufacturer's warranty.

| EXPLICATION DES SYMBOLES : | | | |
|--|--------------------------------------|--|-------------------------|
|  | NON STERILE |  | BATCH NUMBER |
|  | MANUFACTURER | REF | REFERENCE |
|  | DEVICE FOR SINGLE USE - DO NOT REUSE |  | CONSULT THE INSTRUCTION |